

# USER MANUAL



Model No. PM60  
PM60E



## SAVE THESE INSTRUCTIONS

**CAUTION**

Federal (USA) law restricts this device to sale by or on the order of a physician.

**PRECISION**  **MEDICAL**

300 Held Drive  
Northampton, PA 18067

Tel: 610-262-6090  
Fax: 610-262-6080

[www.precisionmedical.com](http://www.precisionmedical.com)

# CONTENTS

RECEIVING/INSPECTION.....	2
INTENDED USE .....	2
READ ALL INSTRUCTIONS BEFORE USING .....	2
EXPLANATION OF ABBREVIATIONS .....	2
SAFETY INFORMATION - WARNINGS AND CAUTIONS ...	2
SPECIFICATIONS .....	4
COMPONENT DESCRIPTION .....	5
OPERATING INSTRUCTIONS.....	6
CLEANING/MAINTENANCE .....	7
RETURNS.....	8
DISPOSAL INSTRUCTIONS .....	8
REPLACEMENT PARTS .....	8
TROUBLESHOOTING.....	9
LIMITED WARRANTY .....	10

## RECEIVING/INSPECTION

Remove the Precision Medical Easy Vac Aspirator from packaging and inspect for damage. If there is any damage, DO NOT USE and contact your Equipment Provider.

## INTENDED USE

To supply a continuous vacuum source adequate within the stated operating vacuum range, to aspirate fluids in medical suctioning procedures by trained healthcare professionals.

## READ ALL INSTRUCTIONS BEFORE USING

This manual instructs a Professional to install and operate the Easy Vac Aspirator. This is provided for your safety and to prevent damage to the Easy Vac Aspirator. If you do not understand this manual, DO NOT USE the Easy Vac Aspirator and contact your Equipment Provider.

## EXPLANATION OF ABBREVIATIONS

inHg	Inches of Mercury
mmHg	Milimeters of Mercury

## SAFETY INFORMATION - WARNINGS AND CAUTIONS

### DANGER

Indicates an imminently hazardous situation which, if not avoided, will result in death or serious injury.

### WARNING

Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

### CAUTION

Indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury.

### CAUTION

Used without the safety alert symbol indicates a potentially hazardous situation which, if not avoided, may result in property damage.



CONSULT ACCOMPANYING DOCUMENTS.



Type B equipment.

## **⚠ DANGER**

### **TO REDUCE THE RISK OF ELECTROCUTION:**

- Always unplug the Easy Vac Aspirator immediately after using.
- **DO NOT** use while bathing.
- **DO NOT** place or store Aspirator where it can fall or be pulled into a tub or sink.
- **DO NOT** place or drop into water or other liquid.
- **DO NOT** reach for a Aspirator that has fallen into water. Unplug immediately.
- **DO NOT** attempt to repair Aspirator. Any attempt could result in electrical shock.
- **DO NOT** use near any type of flame or flammable/explosive substances.
- **DO NOT** allow Aspirator to exceed Operating or Storage Temperature specifications.
- Caution must be exercised when using this Aspirator where oxygen is being administered.

## **⚠ WARNING**

### **TO REDUCE RISK OF BURNS, ELECTROCUTION, FIRE, OR INJURY TO PERSONS:**

- The Easy Vac Aspirator should never be left unattended when running.
- Close supervision is necessary when the Aspirator is used by, on, or near children or invalids.
- Use the Easy Vac Aspirator only for its "Intended Use" as described in this manual.
- **DO NOT** use attachments not recommended by the Precision Medical, Inc.
- **DO NOT** use when sleepy or drowsy.
- **DO NOT** use if Easy Vac Aspirator has a damaged cord or plug.
- **DO NOT** use if Easy Vac Aspirator has been dropped or damaged.
- **DO NOT** cover air vents. Restriction of air flow will cause Aspirator to overheat, causing the thermal reset to shut down the Aspirator.
- **DO NOT** drop or insert any object into any opening or hole.
- Connect this Aspirator to a properly grounded outlet only. See Grounding Instructions.
- Keep the cord away from heated surfaces.

## ⚠CAUTION

If for any reason any label becomes illegible or lost, contact Precision Medical, Inc..

## CAUTION

- Inspect the Easy Vac Aspirator for visual damage, DO NOT USE if damaged.
- Ground reliability can only be achieved when the power cord is connected to an equivalent receptacle marked "Hospital Grade".

## SPECIFICATIONS

**Operating Vacuum Range:** 0-26 inHg (0-660 mmHg)

**Dimensions:**

Length:	13.25 in	(33.7 cm)
Width:	8.00 in	(20.3 cm)
Height:	9.75 in	(24.8 cm)

**Weight:** 12.6 lbs (5.7 kg)

**Shipping Weight:** 14.5 lbs (6.6 kg)

**Electrical:** PM60: 108 -132 VAC, 60 Hz, 5 Amps  
PM60E: 207 - 253 VAC, 50Hz, 2.4 Amps

**Operating Temperature Range:** 50°F to 104°F (10°C to 40°C)

**Storage Temperature Range:** -4°F to 104°F (-20°C to 40°C)

**Filter:** Inline, Hydrophobic

**Safety Systems:** Automatic Resetting Thermal Limit Switch

### Equipment Classification:

Class I	- with respect to protection from electrical shock
Type B	- degree of protection against electrical shock
IPX0	- degree of protection against ingress of liquids
Mode of Operation	- continuous

The Easy Vac Aspirator meets the EMC requirements of EN 550011:1998.


The Easy Vac Aspirator is not suitable for use in the presence of a flammable anesthetic mixture with air, oxygen or nitrous oxide.

Specifications are subject to change without prior notice.


# COMPONENT DESCRIPTION

## ⚠ CAUTION

If for any reason any label becomes illegible or lost, contact Precision Medical, Inc..



Standard 800cc  
Suction Canister



# GROUNDING INSTRUCTIONS

The Easy Vac Aspirator should be grounded. In the event of an electrical short circuit, grounding reduces the risk of electrical shock by providing an escape wire for the electrical current. This Aspirator is equipped with a cord having a grounding wire with a grounding plug. The plug must be plugged into an outlet that is properly installed and grounded.

## **⚠ DANGER**

Improper use of the grounding plug can result in a risk of electrical shock.

## **⚠ WARNING**

Ground reliability can only be achieved when the power cord is connected to an equivalent receptacle marked "Hospital Grade".

# OPERATING INSTRUCTIONS

## **CAUTION**

Inspect the Easy Vac Aspirator for visual damage, **DO NOT USE** if damaged.

1. Place Easy Vac Aspirator near a properly grounded electrical receptacle on a table, floor or other horizontal, smooth, hard surface.  
Allow a minimum of 4in (10 cm) clearance on all sides to ensure adequate air flow for cooling of Aspirator.
2. Press down on Aspirator until Suction Feet are secure.
3. Be sure the On/Off Switch is in the "OFF" position.
4. Be sure the Inline Hydrophobic Filter is properly connected to the vacuum port on the Suction Canister lid.
  - The Inline Hydrophobic Filter prevents liquid from entering the vacuum pump.

## **⚠ CAUTION**

**DO NOT** run the Easy Vac Aspirator without Inline Hydrophobic Filter. This will void the warranty.

5. Plug the Aspirator into a properly grounded receptacle.
6. Turn the On/Off Switch to the "ON" position.

7. Block the patient connection on the Suction Canister and adjust the Vacuum Control Knob to the desired vacuum.  
To INCREASE vacuum - Turn knob CLOCKWISE  
To DECREASE vacuum - Turn knob COUNTERCLOCKWISE
  - Vacuum should be set as prescribed by Physician or Healthcare Professional.
8. Connect one end of the Patient Suction Tubing to the port on the Canister marked Patient Tubing.
9. Connect the other end of the Patient Suction Tubing to the suction catheter and begin procedure as instructed by Healthcare Provider.
10. When procedure is completed, press the On/Off Switch to the "OFF" position.
  - DO NOT overfill the canister
11. Empty the Suction Canister as follows:
  - Remove the canister from its bracket and empty as instructed by your Physician or Healthcare Professional.

### **CAUTION**

Canister is for single patient use only.

12. Replace Suction Canister and reconnect all suctioning tubing.

## **CLEANING/MAINTENANCE**

### **⚠WARNING**

Turn the Easy Go Vac Aspirator off and disconnect from any external power source before cleaning.

### **Inline Hydrophobic Filter:**

This filter cannot be cleaned. If it becomes contaminated or clogged, it must be replaced.

#### **To replace Inline Hydrophobic Filter:**

1. Disconnect filter assembly from Suction Canister and remove Suction Canister from the Aspirator.
2. Remove filter assembly from the Aspirator.
3. Dispose of filter.
4. Install new filter (Part #1690) and reconnect tubing and elbow.

**Tubing** is single patient use only. DO NOT clean.



## Cleaning for Exterior Components:

Use a cloth or sponge dampened with soap and water.  
Wipe dry with a clean cloth or paper towel.

## RETURNS

Returned products require a Returned Goods Authorization (RGA) number, contact Precision Medical, Inc. All returns must be packaged in sealed containers to prevent damage. Precision Medical, Inc. will not be responsible for goods damaged in transit. Refer to Precision Medical, Inc. Return Policy available on the Internet, [www.precisionmedical.com](http://www.precisionmedical.com).


## DISPOSAL INSTRUCTIONS

Dispose of Bio-waste products as described by your Healthcare Professional or Healthcare Facility protocol.

### **⚠WARNING**

If Aspirator is used for treatment of infectious disease, consult your Physician or Healthcare Professional for recommended procedures for proper disposal.

**Please Recycle**



## REPLACEMENT PARTS

<b>Description</b>	<b>Part #</b>
Patient Suction Tubing	1955
Filter	1690
Filter Assembly	502012
Elbow Canister Outlet	1616
800cc Canister, (Case-10 ea.)	502519-10

# TROUBLESHOOTING

If the Easy Vac Aspirator fails to function, consult the Troubleshooting Guide.

If the problem cannot be corrected, consult your Equipment Provider.

<b>Problem</b>	<b>Probable Cause</b>	<b>Remedy</b>
Aspirator will not operate	<ol style="list-style-type: none"> <li>1. On/Off switch in "OFF" position</li> <li>2. Power Cord not connected properly</li> <li>3. Thermal limit exceeded</li> </ol>	<ol style="list-style-type: none"> <li>1. Push On/Off switch to "ON" position</li> <li>2. Check plug at outlet</li> <li>3. Shut Aspirator off, wait 15 minutes, be sure vents are clear</li> </ol>
No vacuum on gauge / no vacuum at Patient Connection Tubing	<ol style="list-style-type: none"> <li>1. Hydrophobic Filter not connected properly</li> <li>2. Suction Canister not connected or assembled properly</li> <li>3. Defective Suction Canister</li> <li>4. Defective Regulator</li> </ol>	<ol style="list-style-type: none"> <li>1. Confirm and tighten connections</li> <li>2. Confirm and tighten</li> <li>3. Replace Suction Canister</li> <li>4. Contact your Equipment Provider</li> </ol>
No vacuum on gauge / vacuum at Patient Connection Tubing	<ol style="list-style-type: none"> <li>1. Defective vacuum gauge</li> </ol>	<ol style="list-style-type: none"> <li>1. Contact your Equipment Provider</li> </ol>
Vacuum on gauge / no vacuum at Patient Suction Tubing	<ol style="list-style-type: none"> <li>1. Clogged Hydrophobic filter</li> <li>2. Canister overfilled</li> </ol>	<ol style="list-style-type: none"> <li>1. Replace Hydrophobic filter</li> <li>2. Empty canister</li> </ol>
Unable to adjust vacuum	<ol style="list-style-type: none"> <li>1. Defective regulator</li> </ol>	<ol style="list-style-type: none"> <li>1. Contact your Equipment Provider</li> </ol>

## **LIMITED WARRANTY AND LIMITATION OF LIABILITY**

Precision Medical, Inc. warrants that EASY VAC Aspirator (the Product) will be free of defects in workmanship and/or material for the following period:

Five (5) years from shipment

Should any failure to conform to this warranty appear within the applicable period, Precision Medical, Inc. shall, upon written notification thereof and substantiation that the goods have been stored, installed, maintained and operated in accordance with Precision Medical, Inc.'s instructions and standard industry practice, and that no modifications, substitutions, or alterations have been made to the goods, correct such defect by suitable repair or replacement at its own expense.

**ORAL STATEMENTS DO NOT CONSTITUTE WARRANTIES.**

The representatives of Precision Medical, Inc. or any retailers are not authorized to make oral warranties about the merchandise described in this contract, and any such statements shall not be relied upon and are not part of the contract for sale. Thus, this writing is a final, complete and exclusive statement of the terms of that contract.

**THIS WARRANTY IS EXCLUSIVE AND IS IN LIEU OF ANY WARRANTY OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE OR OTHER WARRANTY OF QUALITY, WHETHER EXPRESS OR IMPLIED.**

Precision Medical, Inc. shall not under any circumstances be liable for special, incidental or consequential damages including but not limited to lost profits, lost sales, or injury to person or property. Correction of non-conformities as provided above shall constitute fulfillment of all liabilities of Precision Medical, Inc. whether based on contract, negligence, strict tort or otherwise. Precision Medical, Inc. reserves the right to discontinue manufacture of any product or change product materials, designs, or specifications without notice.

Precision Medical, Inc. reserves the right to correct clerical or typographical errors without penalty.

# *PRECISION* *MEDICAL*®

Precision Medical, Inc.  
300 Held Drive  
Northampton, PA 18067 USA

Toll Free Tel: 800-272-7285  
Toll Free Fax: 800-353-1240  
Tel: (+001) 610-262-6090  
Fax: (+001) 610-262-6080

[www.precisionmedical.com](http://www.precisionmedical.com)

**ISO 13485 Certified**