

**NBN INFUSIONS**  
**GENERAL CARE OF THE MONITOR**

<b><u>LEARNING OBJECTIVE</u></b>	<b><u>EXPLANATION</u></b>
Please observe the following guidelines for care and cleaning of the monitor.	<ol style="list-style-type: none"><li>1. <u>Never</u> immerse the monitor cords, or wires in water.</li><li>2. <u>Never</u> use cleaners, alcohol, etc. to clean the monitor. Use only a clean, damp cloth to wipe off the monitor and cords and <u>only when the monitor is unplugged.</u></li><li>3. <u>Never</u> remove any knobs, buttons, or screws from the monitor.</li><li>4. <u>Never</u> attempt any repairs on the monitor, nor allow anyone except authorized NBN Infusion personnel to repair the unit.</li><li>5. <u>Never</u> place the monitor on top or next to any other electrically operated equipment such as clock radios, stereos, TV's, etc.</li><li>6. <u>Never</u> use the monitor while bathing your child.</li><li>7. <u>Never</u> place the monitor in an area where it could be reached by untrained individuals, i.e., children.</li><li>8. <u>Never</u> expose the monitor to extreme temperature, i.e., on top of a radiator, heater, etc.</li><li>9. <u>Never</u> allow the electrical cord to lie in the manner that it may be stepped on or tripped over.</li><li>10. <u>Never</u> tamper with the three-prong plug on the electrical cord. If your outlets at home do not accept a three-prong plug, buy and properly ground a three-prong adaptor.</li><li>11. <u>Always</u> hand-carry the Monitor during bus, taxi, plane, train, and car trips. <u>Never</u> allow the Monitor to be Transported with other "baggage".</li><li>12. The Monitor <u>always</u> must be placed on a sturdy, level surface.</li></ol>