

INOGEN ONE

Why am I smiling?





I weigh less than your cat.

No more worries about running out of oxygen ever again.

No more relying on others for your oxygen supply.

And no more deliveries.

The INOGEN ONE is the first oxygen system designed to provide unparalleled freedom for the active oxygen user. With the INOGEN ONE, you'll no longer need separate pieces of equipment for home, portable, or travel use. Best of all, because it gets oxygen from the air around you, there are no gas or liquid tanks to fill. That means your INOGEN ONE will never run out of oxygen.

The INOGEN ONE is made to make your life easier. It's light enough to carry over your shoulder, and it comes with a cart to roll behind you when you're on the go. We like to think of it as a "polite companion"... and since it's whisper quiet most people don't know it's running.

I'm smaller than your coffee maker.

The INOGEN ONE fits virtually anywhere. Now you can reclaim that room in your home occupied by a noisy concentrator, spare oxygen tanks, or a large liquid oxygen reservoir. The INOGEN ONE can sit discreetly at your feet, behind a chair, next to the bed, or over your shoulder. All of your oxygen needs can be met with this one small device - day or night, home or away.

Technology not only makes the INOGEN ONE smaller and lighter, it makes it easier to operate too. Just a couple of buttons. That's it. There are no regulators to switch or liquid oxygen freeze-ups. Once you try it, you'll wonder why anyone would use something else. Embrace the simplicity of the INOGEN ONE and see how it can make your life easier.



INOGEN ONE size: 11.62" L x 6.00" W x 12.39" H
(29.51cm x 15.24cm x 31.47cm)

www.inogenone.com



I can go wherever you want.

Head home when you *want* to, not when your oxygen gauge says you *have* to.

Restoring your freedom and control of your day is a big part of what makes the INOGEN ONE so special. You can power your INOGEN ONE by plugging it into a wall outlet, DC outlet in your car, or for hours at a time with the included battery. So go ahead, jump in the car to run errands or take a weekend trip to see family. Come and go as you want without having to switch equipment. Or take that long trip you thought you'd never again have the chance to take. The INOGEN ONE is made to follow you anywhere, whether it's around town, on the road, on a cruise, or in the air. Where you go is now up to you!

Speaking of air travel, many airlines now allow you to use the INOGEN ONE when you fly.

Now you can take your oxygen with you, avoiding the hassle of arranging for oxygen at the airport with your provider and the expense of airline-supplied oxygen during flight. Visit www.inogen.net/faa for a list of airlines that currently allow the INOGEN ONE.

oxygen.anytime.anywhere.

What if you could have one simple oxygen system that was versatile enough to use anytime, anywhere? What if you didn't have to wait for deliveries anymore, or spend time arranging for travel oxygen? What if you could simplify your oxygen needs and have one single solution that fit your life?

Getting your INOGEN ONE oxygen system is quick and easy. If you have a home healthcare provider already, call them and ask for one by name. Accept no substitutes.

If this is your first home oxygen device, or if you're having trouble finding the Inogen One, call us and we'll help you find a supplier near you.

Here's a list of accessories included with your INOGEN ONE:



AC Power Supply

When plugged in, the INOGEN ONE delivers an unlimited supply of oxygen.



Mobile Power Charger

Use this to power your INOGEN ONE in your car, boat, or RV.



Lithium Ion Battery

One battery lasts for up to three hours. When it's plugged in, the battery automatically recharges.



Cart & Carry Bag

Designed specifically for the INOGEN ONE, take advantage of your new independence and mobility.

Technical Specifications

Weight 8.31lbs (3.7kg) with Battery: 9.8lbs (4.4kg)

Size Length: 11.62" (29.51cm) Width: 6.00" (15.24cm)
Height including handle: 12.39" (31.47cm)

Oxygen Flow Pulse dose delivery system. Five flow settings from 1-5

Power Sensing AC power supply 100-240V, 50-60Hz (auto sensing to allow worldwide use)
Mobile Power Charger for mobile DC use (will recharge battery)

Battery Duration 3 hours*, recharge approximately 3 hours with AC or DC power

Noise Level 37.4 Decibels*

Warranty 3 year limited warranty**

Operation Simple control functions and easy-to-read LCD display

Use Stationary, portable and travel. Clinically effective for day and night use.

* At setting 2, in typical home use.

** See limited warranty statement for complete details.



Manufactured by
326 Bollay Drive
Goleta, California 93117, USA
805-562-0515
www.inogen.net
info@inogen.net



*For the latest in community and
lifestyle information for people
living with oxygen therapy,
visit www.oxygenation.com.*

Europe Authorized Representative
HealthLink Europe BV
Centaurusweg 123
NL-5015 TC Tilburg
The Netherlands
Tel: +31 (0)13 547 9364
Fax: +31 (0)13 547 9302